

LF-1, Drawer unit

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS

Drawer internal dimensions (W×D×H)	235 × 385 × 50 mm</br>235 × 385 × 150 mm
Overall dimensions (W×D×H)	300 × 450 × 705 mm
Weight	28 kg




CAT. NUMBER

BS-050101-BK

LF-1, Drawer unit

ACCESSORIES



UVC/T-AR

BS-040102-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVC/T-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, plexiglas walls and ...

[read more](#)



UVC/T-M-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVC/T-M-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)




LT-120

BS-040107-EK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



UVT-B-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVT-B-AR is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



UVT-S-AR

BS-040107-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVT-S-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



LT-150

BS-040107-JK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



LT-180

BS-040107-HK

Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



T-4

BS-040101-BK

Table

New modular design of laboratory furniture provides flexibility and ease of use.



T-4L

BS-040107-BK

Table

New modular design of laboratory furniture provides flexibility and ease of use.

LT-120, Height adjustable desk

DESCRIPTION

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety of laboratory settings. The worktop is easy to clean and maintain, ensuring a hygienic work environment. The table is also adjustable in height, which allows you to customize the table to your preferred working height for added comfort and convenience. In addition to its standard features, this laboratory table also offers an optional drawer with an anti-slippery surface and a retractable block of electrical sockets. This added feature allows for convenient storage of lab equipment and tools, and the anti-slippery surface ensures that items stay securely in place. The retractable block of electrical sockets allows for easy access to power outlets, eliminating the need for extension cords and reducing clutter on the work surface. This table is an all-in-one solution for any laboratory or research facility, providing a sturdy, reliable and practical work surface that can stand up to the demands of daily use.

- Standard delivery set includes 2 legs with 2 motors.
- Controller with memory - 4 memory modes, digital height indicator, button with up/down movement.
- Adjustable height: 625 - 1280 mm.
- Lifting capacity: 89 kg.
- Work surface from chemically resistant HPL, width - 1200mm, depth-600mm, thickness of 25 mm, blue edge.

Drawer and retractable socket block comes as an accessories.



CAT. NUMBER

BS-040107-EK

LT-120, Height adjustable desk

SPECIFICATIONS

Maximum load	89 kg
Overall dimensions (W×D×H)	1200 × 600 × 625 - 1280 mm
Weight	45 kg

ACCESSORIES



LTD-1
 BS-040107-OK
 Drawer

Drawer with non slip silicone mat and soft close function. Can be mounted afterwards.



SBR-1
 BS-040107-NK
 Socket block

Socket block, retractable, with blue LED.



LF-1
 BS-050101-BK
 Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.

LT-150, Height adjustable desk

DESCRIPTION

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety of laboratory settings. The worktop is easy to clean and maintain, ensuring a hygienic work environment. The table is also adjustable in height, which allows you to customize the table to your preferred working height for added comfort and convenience. In addition to its standard features, this laboratory table also offers an optional drawer with an anti-slippery surface and a retractable block of electrical sockets. This added feature allows for convenient storage of lab equipment and tools, and the anti-slippery surface ensures that items stay securely in place. The retractable block of electrical sockets allows for easy access to power outlets, eliminating the need for extension cords and reducing clutter on the work surface. This table is an all-in-one solution for any laboratory or research facility, providing a sturdy, reliable and practical work surface that can stand up to the demands of daily use.

- Standard delivery set includes 2 legs with 2 motors.
- Controller with memory - 4 memory modes, digital height indicator, button with up/down movement.
- Adjustable height: 625 - 1280 mm.
- Lifting capacity: 86 kg.
- Work surface from chemically resistant HPL, width - 1500mm, depth-600mm, thickness of 25 mm, blue edge.

Drawer and retractable socket block comes as an accessories.

SPECIFICATIONS

Maximum load	86 kg
Overall dimensions (W×D×H)	1500 × 600 × 625 - 1280 mm
Weight	49 kg

ACCESSORIES



LTD-1
 BS-040107-OK
 Drawer

Drawer with non slip silicone mat and soft close function. Can be mounted afterwards.



SBR-1
 BS-040107-NK
 Socket block

Socket block, retractable, with blue LED.



LF-1
 BS-050101-BK
 Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.



CAT. NUMBER

BS-040107-JK

LT-150, Height adjustable desk

LT-180, Height adjustable desk

DESCRIPTION

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety of laboratory settings. The worktop is easy to clean and maintain, ensuring a hygienic work environment. The table is also adjustable in height, which allows you to customize the table to your preferred working height for added comfort and convenience. In addition to its standard features, this laboratory table also offers an optional drawer with an anti-slippery surface and a retractable block of electrical sockets. This added feature allows for convenient storage of lab equipment and tools, and the anti-slippery surface ensures that items stay securely in place. The retractable block of electrical sockets allows for easy access to power outlets, eliminating the need for extension cords and reducing clutter on the work surface. This table is an all-in-one solution for any laboratory or research facility, providing a sturdy, reliable and practical work surface that can stand up to the demands of daily use.

- Standard delivery set includes 2 legs with 2 motors.
- Controller with memory - 4 memory modes, digital height indicator, button with up/down movement.
- Adjustable height: 625 - 1280 mm.
- Lifting capacity: 83 kg.
- Work surface from chemically resistant HPL, width - 1800mm, depth-600mm, thickness of 25 mm, blue edge.

Drawer and retractable socket block comes as an accessories.



SPECIFICATIONS

Maximum load	83 kg
Overall dimensions (W×D×H)	1800 × 600 × 625 - 1280 mm
Weight	52 kg

ACCESSORIES



LTD-1
 BS-040107-OK
 Drawer

Drawer with non slip silicone mat and soft close function. Can be mounted afterwards.



SBR-1
 BS-040107-NK
 Socket block

Socket block, retractable, with blue LED.



LF-1
 BS-050101-BK
 Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.

CAT. NUMBER

BS-040107-HK

LT-180, Height adjustable desk

T-4, Table

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS


Maximum load	50 kg
Overall dimensions (W×D×H)	800 × 600 × 770 mm
Weight	23 kg



CAT. NUMBER

BS-040101-BK	T-4, Table
--------------	------------

ACCESSORIES



UVC/T-AR

BS-040102-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVC/T-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, plexiglas walls and ...

[read more](#)




UVC/T-M-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVC/T-M-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



UVT-B-AR

With built in socket and inlet for power cords
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box
UVT-B-AR is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



LF-1

BS-050101-BK
Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.

T-4L, Table

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS

Maximum load	75 kg
Overall dimensions (W×D×H)	1290 × 600 × 770 mm
Weight	36 kg



CAT. NUMBER

BS-040107-BK

T-4L, Table

ACCESSORIES



UVT-S-AR
BS-040107-AAA
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-S-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



LF-1
BS-050101-BK
Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.